

INTENDED USE

The VIRGO® FTA-ABS IgG indirect fluorescent antibody (IFA) test is intended for the detection and titration of *Treponema pallidum* IgG antibodies in human sera. The test procedure is based on the method developed by Hunter, *et al.*, to remove nonspecific reactions by absorbing test samples with nonpathogenic treponemes. The current procedure is described in the 1990 Manual of Tests for Syphilis.¹

SUMMARY

Treponema pallidum, the etiological agent of syphilis, produces in the host two types of antibodies: 1) nontreponemal antibodies (regain), which react with lipid antigens and 2) treponemal antibodies, which react with *T. pallidum* and closely related strains.²

Detection of regain can be successfully accomplished by use of various commercially available tests. Unfortunately, acute or chronic infections such as malaria, leprosy, infectious mononucleosis and upper respiratory diseases as well as collagen and immunological diseases such as rheumatoid arthritis and lupus erythematosus can produce false positive regain tests.^{1,3-12}

Tests for syphilis employing treponemal antigens are of most value in testing sera from patients presenting diagnostic problems. Such individuals most frequently have reactive regain tests in the absence of clinical or historical evidence of syphilis, or may have nonreactive regain tests and clinical signs of late syphilis. A reactive treponemal test is considered good evidence of past or present syphilitic infection, providing that other treponematoses can be ruled out.


The FTA-ABS kit manufactured by Hemagen Diagnostics, Inc., provides all the necessary reagents for the rapid determination of antibody to *Treponema pallidum* in human sera. Antigenic substrate, control sera, FITC conjugate, sorbent, buffer, coverslip mounting media and an instructional insert are included in the kit.

PRINCIPLE OF THE TEST


The VIRGO fluorescent antibody assays utilize the indirect method of fluorescent antibody staining, first described by Weller and Coons in 1954.¹⁶ The procedure is carried out in two basic reaction steps. In step one, the human serum to be tested is brought into contact with the antigenic substrate. Antibody, if present in the test serum, will attach to the antigen, forming an antigen-antibody complex. If the serum being tested does not contain antibody for this particular antigen, no complex is formed and all the serum components are washed away in the rinse step. The second step involves adding a fluorescein labeled anti-human antibody to the test wells. If the specific antigen-antibody complex is formed in step one, the fluorescein labeled antibody will attach to the antibody moiety of the complex in step two. A positive reaction, bright apple-green fluorescence, can be seen with the aid of a fluorescence microscope.

Principle of Indirect Fluorescent Antibody Testing

Step 1



Step 2



CONTENTS OF THE KIT

902300	902301	Test Kit Product Number
60 Tests	200 Tests	Number of Tests per Kit
15 x 4 well	25 x 8 well	Slides: Fixed <i>Treponema pallidum</i> (Nichols Strain)
1 (10 mL)	2 (10 mL)	Vial(s) Sorbent (liquid): A concentrated, standardized extract of a Reiter Treponeme culture. Sorbent removes the nonspecific antibody for treponemes in

1	1	human serum without significantly reducing the specific antibodies to <i>T. pallidum</i> .
1	1	Vial FTA-ABS Reactive Control: Lyophilized human serum
1	2	Vial FTA-ABS Nonspecific Control: Lyophilized human serum
1	2	Vial(s) FITC Conjugate: Liquid, inactivated rat anti-human IgG (heavy & light chains)
2	4	Vial(s) FTA-ABS Conjugate Diluent contains 2% Tween-80
1	1	*Packages Powdered Phosphate Buffer: (PBS) pH 7.4 ± 0.2
2	3	*Vial(s) (2 mL) Buffered Glycerol
		*Packages (5 each) Blotters

*These components may be interchanged between different master lots. Additional supplies are available from Hemagen Diagnostics, Inc.

MATERIALS REQUIRED BUT NOT SUPPLIED

Test tubes and racks for making dilutions

Pipettes for preparing dilutions

Coverslips, 22 x 50 mm., No. 1 thickness

Humidified chamber

Magnetic stir plate (optional)

Staining dish and slide-holder rack

Fluorescence microscope. Refer to manufacturer's instruction manual for the filter system that gives optimum results for FITC (Maximum excitation wavelength = 490 nm. Mean emission wavelength = 520 nm.)

PRECAUTIONS

1. HANDLE ALL ASSAY SPECIMENS, SLIDES, REACTIVE AND NON-SPECIFIC CONTROLS AS IF CAPABLE OF TRANSMITTING INFECTIOUS AGENTS.

All human blood components of the kit have been tested by approved FDA methods and found to be negative for both hepatitis B surface antigen (HBsAg) and for antibodies to human immunodeficiency virus type 1. Because no test method can offer complete assurance that HIV, hepatitis B virus, or other infectious agents are absent, specimens and kit reagents should be handled at the Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen.^{17,18}

- The antigenic substrates are fixed in methanol.
- Do not pipette by mouth.
- Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
- All materials used in this assay, including reagents, samples and wiping materials should be disposed of in a manner that will inactivate infectious agents.

HANDLING PRECAUTIONS

- For In Vitro Diagnostic Use.** No U.S. Standard or Potency.
- Do not use the kit or individual reagents beyond their labeled expiration dates.
- The components of this kit have been tested as a unit. **Do not** interchange components from other sources or from different master lots, except as noted.
- Protect the conjugate from prolonged exposure to light.

REAGENT STORAGE AND STABILITY

- T. pallidum* Antigen Slides** – store at 2-8°C
- Reactive Control, Nonspecific Control, Sorbent and Conjugate** – store lyophilized reagents and sorbent at 2-8°C. **After rehydration**, Reactive Control, Nonspecific Control and FITC Conjugate should be made up in aliquots and stored at -20°C or colder if not used within one week.
- Powdered PBS, Diluent and Buffered Glycerol** can be stored at 2-30°C if desired. The test kit can be used through the expiration date on the outer box label.

NOTE: Precautions were taken in the manufacture of this product to protect the reagents from contamination. After reconstitution, care should be exercised to protect the reagents in this kit from contamination. If constant storage temperature is maintained, reagents and substrate will be stable for the dating period of the kit.

REAGENT PREPARATION

Allow reagents and slides to reach room temperature 15 to 30 minutes before use.

- Slides and Glycerol:** Ready to use.
- PBS:** Dissolve contents of one package in 1 liter of distilled or deionized water. Seal container to prevent contamination or evaporation.
- Controls:** Dehydrate the Reactive and Nonspecific Controls with 1.0 mL of distilled or deionized water. Aliquot for storage at -20°C if not used within one week.
- Conjugate:** Dehydrate with 2.0 mL of FTA-ABS Conjugate Diluent (supplied). Aliquot for storage at -20°C or colder if not used within one week.

SPECIMEN COLLECTION AND HANDLING

- Serum samples may be stored at room temperature for up to 24 hours. For longer term storage, they may be stored at 2-8°C (for up to three days), or frozen at -20°C or colder. Place at 37°C only until the samples are thawed. Remove and mix thoroughly before use. Self-defrosting freezers are not recommended. Avoid multiple free-thaw cycles.¹⁹
- Optimal performance of the VIRGO FTA-ABS IFA depends upon the use of fresh serum samples. Specimens should be collected aseptically. Early separation from the clot prevents hemolysis of serum.²⁰ No anticoagulants or preservatives should be added.
- For best results, another sample should be drawn if bacteriological contamination or lipids are present. If another sample cannot be obtained, filtering (0.45µ) or centrifugation (approximately 3000 x G for 10 minutes) is required.
- Occasionally, the specimen may contain certain proteolytic enzymes which attack and digest the substrate. This is especially true of specimens contaminated with microorganisms. Such specimens may be heated to 56°C for 30 minutes. If this fails to reduce the enzymatic activity, another sample should be obtained from the patient.

TEST PROCEDURE

Specimens may contain infectious agents and should be handled accordingly.

For optimal results, DO NOT allow substrate wells to dry out while performing the test.

- Remove the slides and required reagents from the refrigerator** and allow them to reach room temperature (15 to 30 minutes).
- Heat all sera (control serum aliquots and test samples) for 30 minutes at 56°C in a water bath prior to testing.** Sera previously heated and stored in aliquots should be heated to 56°C for at least 10 minutes before reuse.¹² Allow the heated sera to cool to room temperature before proceeding to next step.
- Seven controls will be required for each group of samples being tested.** The following table indicates the dilution for each control:

Controls	Dilution
1. Reactive Control	1:5 in PBS
2. Reactive Control	1: 5 in Sorbent
3. Nonspecific Control	1:5 in PBS
4. Nonspecific Control	1:5 in Sorbent
5. Minimally Reactive Control to Yield a 1+ Fluorescent Reaction	See FTA-ABS Dilution Notice
6. Nonspecific Fluorescence Controls	Undiluted PBS
7. Nonspecific Fluorescence Controls	Undiluted Sorbent

These dilutions are prepared by placing 0.1 mL of the diluent in a test tube and adding 0.025 mL of the control to give a 1:5 dilution. Dilutions must be thoroughly mixed. (See Quality Control Section).

- Preparation of 1:5 dilutions of test samples:** To appropriately marked tubes, add 0.1 mL of sorbent. Next add 0.025 mL of previously heated (but now cooled) sample. Mix by pipetting. Incubate at room temperature for 15 minutes. During incubation, mix by shaking the test tube rack a few times.
- Remove the slides from the pouch** just before use and label, keeping in mind that seven controls must be run each time the test is performed.
- Cover each well with one sample, either a control serum dilution or a dilution of a patient serum specimen (see Step #3 above).** (~10 – 20 µL per well).
- Incubate in a humidified chamber at 37 ± 2°C for 30 minutes.**
- Rinse the slides briefly** in a light stream of PBS. Do not direct the stream into the wells.
- Rinse the slides thoroughly for 7 minutes** in a staining dish of PBS. **Change the buffer and wash for an additional 8 minutes.** Handle slides gently. Gentle agitation of the buffer is necessary for efficient slide washing.
- Blot the painted mask of the slide** with the blotters provided. **Do not allow the wells to dry before conjugate addition.**
- Cover each well with one drop (~10 µL) of FITC Conjugate.**
- Incubate in a humidified chamber at 37 ± 2°C for 30 minutes.** Protect from intense light.
- Repeat steps 8 and 9. Blot the painted mask of the slide** with the blotters provided. **Do not allow the wells to dry before the addition of conjugate.**
- Place a small drop of Buffered Glycerol in each well** and cover with a coverslip.
- For best results, **the slides should be read immediately at a magnification of 200-500X.** Alternatively, the slides may be read within 24 hours. However, they should be stored at 2-8°C in the dark, and sealed to prevent the mounting fluid from drying.

CRITERIA FOR GRADING FLUORESCENCE INTENSITY

Reading	Intensity of Fluorescence
2+ to 4+	Moderate to Strong
1+	Equivalent to Minimally Reactive (1+) control
± to < 1+	Visible staining but less than (1+)
-	None or vaguely visible but without distinct fluorescence
Atypical	Varied: treponemes appear to be "moth eaten" or to have "beads" of fluorescence throughout their length

INTERPRETATION OF SAMPLE RESULTS

Reporting System for FTA-ABS Test		
Initial Test Reading	Repeat Test Reading	Report
4+		Reactive
3+		Reactive
2+		Reactive
1+ ^a	>1+	Reactive
1+	1+	Reactive Minimal ^b
1+	<1+	Nonreactive
<1+		Nonreactive
N		Nonreactive
Bead fluorescence		Atypical fluorescence observed ^c

^aRetest all specimens that initially give 1+ intensity of fluorescence.

32. **Hunter, EF, MR Adams, LH Orrison, BJ Pender and SA Larsen.** 1979. Problems affecting performance of the fluorescent treponemal antibody-absorption test for syphilis. *J. Clin. Microbiol.*, **9**:163-166.
33. **Julian, AJ, LC Logan and LC Norins.** 1969. Early syphilis: immunoglobulins reactive in immunofluorescence and other serologic tests. *J. Immunol.* **102**:1250.
34. **Hunter, EF, JF Smith and JR Lewis.** 1972. Immunochemical parameters of commercial conjugates for the fluorescent treponemal antibody-absorption test. *Infect. and Immun.*, **5**:858-862.
35. **Wilkinson, AE and P Rodin.** 1976. IgM-FTA test in syphilis in adults in relation to clinical finding. *Brit. J. Vener. Dis.*, **52**:219.
36. **Hardy, PH, Jr. and EE Nell.** 1972. A study of properties of sorbent, the reagent employed in the fluorescent treponemal antibody-absorption (FTA-ABS) test. *Health Lab. Sci.*, **4**:5.
37. **Olitzky, I and RG Blaker.** 1974. Laboratory aids for the diagnosis and treatment of gonorrhoea and syphilis. *Bioscience Laboratory Manual*, **9**.
38. **Deacon, WE, VH Falcone and A Harris.** 1957. A Fluorescent test for treponemal antibodies. *Proc. Soc. Exp. Biol. Med.*, **96**:477-480.
39. **Hunter, EF, WE Deacon and PE Meyer.** 1964. An improved FTA test for syphilis, the absorption procedure (FTA-ABS). *Pub. Health Rep.*, **79**:410-412.
40. **Stout, GW, DS Kellogg, Jr., VH Falcone, BE McGrew and JS Lewis.** 1967. Preparation and Standardization of the sorbent used in the fluorescent treponemal antibody-absorption (FTA-ABS) test. *Health Lab. Sci.*, **4**:5.

SUMMARY OF VIRGO FTA-ABS

IMPORTANT: It is recommended that one be familiar with the detailed procedure in the package insert before using this summary.

Heat inactive serum samples by heating at 56°C for 30 minutes. Allow same to cool.



Dilute samples and control with sorbent or PBS.



Incubate at 23 ± 2°C for 15 minutes while shaking.



Cover wells with patient samples or controls as described.



Incubate slides in a humidified chamber at 37 ± 2°C for 30 minutes.



Rinse slides briefly with PBS. Wash for 15 minutes, changing buffer 1X. Blot slides.



Cover each well with FITC Conjugate.



Incubate slides in a humidified chamber at 37 ± 2°C for 30 minutes.



Repeat wash step described above.



Place a small drop of Buffered Glycerol on each well and cover with a coverslip.



Read slides immediately at 200-500X magnification on fluorescence microscope.

Fluorescent Treponemal Antibody-Absorption/ FTA-ABS IgG IFA

Immunofluorescence Test Kit for the Detection of Antibody to *Treponema pallidum* In Human Sera

FOR *IN VITRO* DIAGNOSTIC USE

Hemagen Diagnostics, Inc. VIRGO® Products Division

Columbia, Maryland 21045

Phone: (800) 436-2436

(443) 367-5500

Web Site: www.hemagen.com



EC Representative
Obelis S.A.
Avenue de Tervueren, 34 Bte 44
B-1040 Brussels, Belgium
tel: (32) 2.732.59.54
fax: (32) 2.732.60.03

P/N 890100-7 Rev. J
November 2006